Colon Transit (Kolonpassagezeit)

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Fabricant: P. & A. Mauch, CH-4142 Münchenstein
CE 0044

Rev. 06/2022

Colon Transit (Colonic Transit Time)

Composition:

Six cellulose capsules, each containing 10 radio-opaque poly- urethane markers containing 40% barium sulphate.

The capsules are numbered and contain markers of the following shapes:

Capsule No./Capsule content:	Capsule No./Capsule content:
Capsule 1: 10 rods	Capsule 4: 10 cubes
Capsule 2: 10 spheres	Capsule 5: 10 small rings
Capsule 3: 10 large rings	Capsule 6: 10 rods

Description:

This test is a diagnostic aid which has been developed in order to investigate the passage of the contents of the colon and the rhythm of elimination. It comprises 6 capsules, which should be swallowed in the numbered sequence and at the given times. Each capsule contains 10 radio-opaque markers which do not undergo any changes in the gastro-intestinal tract. Without any special preparation of the patient, by means of a single X-ray picture (100KV) after 7 days (144 hrs) it is possible to reconstruct the activity of the colon during this period. The calculation is made according to Metcalf et al. (Gastroenterology 1987, 92 40-47) and Chaussade (Gastronétrol. clin. biol. 1986, 10: 385-389):

COLONIC TRANSIT TIME = (Sum of the markers x [time between ingestions / number of the markers per capsule]): **CTT (hours) = sum of markers x 2.4**

Determination of the colonic transit time CTT:

Number of markers in the colon	<=10	11-20	21-30	31-40	41-50	51-
Corresponds to passage time in hours	<=24	26-48	50-72	74-96	98-120	122-
Normal values WOMEN (hours)			<=70			
Normal values MEN (hours)			<=60			

Ingestion:

- Days 1-6: 1 capsule at 9.00 a.m., with breakfast
Day 7: Plain abdominal radiograph at 9.00 a.m.

Instructions for use:

Prescription and explanation of use by the attending specialist. There must be a justifiable indication for taking an X-ray. Safety and efficacy in children (<18 years) have not been established. After oral administration of barium sulphate (B.) in quantities of up to 100 g, side effects have been described only very rarely despite many years of use. These and also the even more rarely reported barium-specific toxic effects could in no case be clearly assigned to a direct effect of pure B. Specific to this product (max. 1.3 g B. in total), no side effects have been reported in this respect in over 30 years of use. If, contrary to expectations, you should experience side effects, contact the attending specialist immediately.




Contraindications:

In case of known hypersensitivity to hydroxypropylmethylcellulose (HPMC / E464), Indigotine - FD&C Blue2 (E132), titanium dioxide (E171), Iron oxide (E172), Barium sulphate (BaSO₄), Pellethane® 2363 TPU.

Warnings:

Not for use in pregnant women due to radiation concerns.

Symbols used:

	Manufacturer		Temperature limit
	Date of manufacture		Follow instructions for use
	Reference number		EU authorised representative
	Batch number		CE marking
	Usable until		Medical device
	Do not use if packaging is damaged		Store in a dry place
	Protect from sunlight		

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